



Highly improved multi-photon imaging with easy-to-use femtosecond laser technology


Image quality in multi-photon microscopy depends drastically on the compression of the laser pulses at the sample. Every optical component in a microscope introduces temporal dispersion that increases the laser pulse duration.

Ensuring optimum pulse compression at the sample is not an easy task, and is harder to do for shorter pulses. Sphere's d-micro was especially developed to solve this problem by simultaneously compressing the pulse and ensuring optimum duration directly at the sample plane, even for the most demanding laser pulses.

The d-micro comprises two COMPACT modules:


- the COMPRESSOR, which allows pre-compensating the dispersion introduced by the microscope
- the MEASURING HEAD, which measures the pulse directly at the sample position

Sphere Ultrafast Photonics can install the d-micro system on a microscope already equipped with a femtosecond laser, or can provide a complete ultra-broadband femtosecond laser with bundled d-micro system, with guaranteed pulse duration below 10 fs at the sample plane.



The images show a greatly improved signal intensity (improved S:N and image contrast) for the same average power level, as well as details that can be seen using a compressed sub 10 fs laser pulse instead of a significantly longer pulse duration.








Credits: Kurt Anderson, Francis Crick Institute / Thorlabs / Sphere

SyncRGB:FLIM

5 um



https://doi.org/10.1364/BOE.10.001891