LED Illumination Covered







- 8 individually controllable channels
- <7 µs TTL triggering
- Software, digital & analogue control



Research Fluorescence Illumination pE-4000

- 16 selectable LED wavelengths
- No modularity
- Compatible with all filter sets



Screening Fluorescence Illumination pE-100

- Instant on/off
- 0-100 % intensity control
- Select from full range of LED wavelengths



Transmitted Applications pT-100

 Optimised transmitted light LED Illumination

LED Illumination Covered





Everyday Fluorescence pE-300 Series

- Broad spectrum
- · Direct or light guide delivery options
- Sustainable ACT label certification

pE-300lite

Simple white light

pE-300^{white}

- Everyday controllable illumination
- USB and TTL control

pE-300^{ultra}

- Fast controllable illumination
- Removable excitation filter holders
- Sequence Runner

Fura-2 Calcium Imaging pE-340^{fura}

- High-speed µs switching
- Enhanced contrast and cell health
- Stable and precise LED illumination

Scan here to find us on WeChat



For more information on how CoolLED products can help you, contact us now:

- t: +44 (0)1264 323040 (Worldwide)
- t: **1-800-877-0128** (USA/Canada) w: **www.CoolLED.com**
- e: info@CoolLED.com

